


BRS System

RADIOLOGICAL SYSTEM WITH FIXED FFD



Multi-purpose system

User friendly

Mobile table

Compact system

Fast and Effective Examinations

BRS SYSTEM

radiological system with fixed FFD

The equipment WHIS RAD is perfect for your routine needs and more. The multi-purpose design allows for a combination of procedures normally performed by much larger and more expensive systems. Easy to use, the swivel arm moves vertically on a column to facilitate upright patient procedures in addition to routine examinations utilising a mobile table.

The B.R.S. enables fast and effective examinations and procedures. Plus it is compact, economic and easy to use.



Technical specifications

Bras


Fixed Source Image Distance	140 cm
Vertical movement	142 cm
Arm rotation	150°

Potter

Bucky's centre/ floor distance	from 45 to 187 cm arm in horizontal position
Skin/ Film distance	< 5 cm
Grid	ratio 12 : 1, 40 lines/ cm, F=150 cm
Unit's dimensions and weight	135 x 159 x 224 cm - 340 kg

Table

Dimensions and weight	Casters table with brakes, DIN Rails for accessories (L x W x H) 200 x 65 x 70 cm - 46 kg
-----------------------	---



Layout example



C/ Pelaya, 13 - Pol. Ind. "Río de Janeiro"
28110 - Algete (Madrid-SPAIN)

Central offices:
Tel: (34) 91 - 628 0544
Fax: (34) 902 190 385
Email: info@sedecal.com

