

PHOENIX



Improvements

- Maximum column travel
- Maximum arm travel
- Status LED indicators

Safety Features

- Proximity sensors
- Bumper



New TFT 9" includes exposure parameters and Integration with the High Frequency X-Ray Generator.



SEDECAL

Calle Pelaya 13 - Pol. Ind. Río de Janeiro
28110 Algete (Madrid - SPAIN)
Central Offices:
Tel: (+34) 91 628 05 44
Fax: (+34) 902 190 385
Email: info@sedecal.com

SEDECAL USA Inc

3190 N. Kennicott Ave.
Arlington Heights, IL 60004
USA
Tel: 1 (847) 394 - 6960
Fax: 1 (847) 394 - 6966
Email: sales@sedecalusa.com



YOUR OEM PARTNER

PHOENIX

you won't believe what you see



Compact, Telescopic, Easy Moving and the most **Powerful** high frequency generator in the market,



YOUR OEM PARTNER



PHOENIX

The new mobile **Phoenix** with telescopic column provides all the technologies attributes used for SEDECAL mobiles with new innovations thinking in the patient and the professional: power, reliability, high quality imaging in digital X-ray, coverage and now also extreme visibility.

The new compact size, the integration of the system plus all the elements added in terms of usability and safety brings Phoenix as the the most complete X-ray mobile in the market.



Fully Digital Integrated

PhoeniX mobile is the new digital generation mobiles from SEDECAL. Includes double TFT totally integrated with the high frequency generator, where the operator will be able to proceed with the exposure, x-ray image preview and image enhancement tools. This mobile has been designed with high advantages compare with the rest of systems nowadays:

- ✓ **Compact Size**
- ✓ **Telescopic Column**
- ✓ **Powerful (Up to 50KW)**
- ✓ **Fully motorized**
- ✓ **Easy Moving**




From excellent to perfect

THE MOST POWERFUL IN THE INDUSTRY

Max.
50 kW

Phoenix is the most powerful high frequency generator mobile in the market, with a range from 20 Kw to 50 Kw Nominal power. The proven reliability of these generators assures **constant performing during the whole day**.

COMPACT SIZE AND TELESCOPIC COLUMN DESIGNED FOR AN ADVANCED VISIBILITY



The new generation of Sedecal mobiles with telescopic column solves any kind of obstacle that users and professional could find while they are driving or parking into the clinic. The new unit provides **maximised vision thanks to the telescopic column and the ergonomic design**. Nurses and doctors will be able to access any location at the hospital in a comfortable and safety way.

MAXIMUM MOBILITY AND COVERAGE



This mobile has been designed with a **brilliant compact size**, which improves maneuverability and access points in any clinic room. **Thanks to the technologies "Easy driven, easy rotation and easy extension"**, the level of mobility and coverage of patient is maximum.

HIGHLY CONFIGURABLE FOR TOP FLEXIBILITY



Phoenix mobile **can be integrated with any wireless detector** and with its digital acquisition software. Provided with new **19" Touch Screen** and multitouch option. **Charging detector batteries are included directly in the mobile**. Also the system offers multiple possibilities to connect external devices for diagnostic