

Ozonette

WOMAN CARE



CE
1639

MD



Technical features

Operating modes syringe automatic loading, continuous mode with independent flow adjustment, Insufflation mode (full user control).

Concentration range 1 to 80 ug/Nml in steps of 1 ug / Nml.

ACM concentration measurement algorithm mathematical calculation.

N Standardization of the measurement (regarding pressure and temperature).

DC dynamic control The system is continuously monitoring the gas flow in order to maintain the concentration.

OC automatic overpressure detection.

Output 3 g/h.

Optimized oxygen consumption. Design oriented to avoid medical oxygen waste.

Compact, high-precision generator tube, inert and highly mechanically strength made of the latest generation of materials (steel 316L).v. Uniform flow distribution.

High voltage generator™.

Ozone removal by double circuit ozone flow Produced, digitally adjustable from 10 to 50 L/h.

Operational control digital touch control.

Graphic LCD screen.

Interface (universal and intuitive).

Weight/Dimensions 3 kg (6.6 lb) / 25 x 27 x 12cm (9.8 x 10.6 x 4.7 inch).

Certificates, CE (IIB class), RoSH.

Compatibility with oxygen cylinder and centralized oxygen supply system.

Universal external power supply medical standard (100 to v240 V / 50 to 60 Hz).

Specific insufflation mode with variable flow adjustment for greater patient comfort and safety with overpressure cut-off under risk analysis.

Ozonette is a compact medical ozone generator that incorporates the latest technological advances. Designed for ozone therapy applications where concentrations must be measured accurately and safely, complying with international standards.

Three modes are available in **Ozonette**:

- Syringe operation is the most advanced on the market with automatic syringe size detection. The ozone output is done by a double valve system and automatic pressure compensation to avoid unwanted leaks ensuring the user and patient safety.
- Continuous ozone mode operation with variable flow setting and overpressure detection is ideal for local bag therapies in both open and closed circuit, or liquid to be ozonized.
- Insufflation mode operation allows to program concentration and variable flow, activating the ozone production and output controlled by the user, providing the volume data and doses supplied during production.

The generator tube is made of high-quality materials, inert to ozone's action and with great mechanical resistance. The high-voltage transformer, specifically designed for our products, ensures optimal ozone generation™ through the dielectric discharge process.vv



SEDECAL S.A.
(Fábrica/Oficinas Centrales)
C/ Pelaya, 9
28110 Algete,
Madrid (SPAIN)
+34 916 280 544
O3@sedecal.com
www.o3.sedecal.com

