



# **PEGASUS**

## **Battery Portable X-Ray System** ✕

Turnkey cost-effective OEM Solution.

# PEGASUS LIGHTWEIGHT OEM SOLUTION

- Battery-powered high-frequency X-Ray generator.
- Robust & lightweight.
- Clinical use.
- Weight < 90 kg..

## TWO VERSIONS

### ANALOG

**READY FOR FUTURE UPGRADE TO DIGITAL.**

- Cassette holder / detector holder.
- Tray for portable PC.




### DIGITAL

- Storage for 35x43 cm or 43x43 cm detector.
- Detector locking. Anti-theft system (optional).
- Detector charging while plugged (optional).



## MOVEMENTS

- Compact design.
- Excellent maneuverability and visibility while driving.
- Ergonomic and easy use.
- Omni-directional wheels.
- Focus-floor distance: 2000 mm
- Column rotation.



Bringing your OEM solution with the greatest patient coverage.  
Pegasus is a real compact, robust & lightweight equipment.

# UP TO 8 KW


## BATTERY DEVICE

# UP TO 20 EXPOSURES



## USER INTERFACE

- Corded handswitch: Preparation, Exposure and collimator light.
- Wireless remote control for Preparation, exposure and collimator light (optional).



# CREATING YOUR X-RAY VISION

Your Strategic OEM Partner 



There has never been an easier way to make a New Product Introduction (NPI) and globally distribute your Product. Find your state-of-the-art solution from the comprehensive Sedecal portfolio.

- Complete product process, from design to transfer-to-production, verification & validation, registration & certification and logistic tracking.
- Highest international quality standards, including Quality management system, CE MARK and FDA CLEARANCE, EMC & electrical safety, X-Ray protection, biocompatibility and sustainability.
- Worldwide registration support, dealing with International organisms and local authorities.
- Heavy investment in R&D, with a cutting-edge engineering and regulatory team, trendy industrial designers and 3D-printing prototyping capability.

Access to all our services and X-ray solutions and SIMPLIFY your NPI on [www.sedecal.com](http://www.sedecal.com)



## SEDECAL S.A.

C/ Pelaya, 9-13  
28110 Algete,  
Madrid, SPAIN

+34 916280544  
[info@sedecal.com](mailto:info@sedecal.com)  
[www.sedecal.com](http://www.sedecal.com)

## SEDECAL USA, INC

7555 N. Caldwell Ave.  
Niles, IL 60714, USA

+001 8473946960  
[salesoem@sedecalusa.com](mailto:salesoem@sedecalusa.com)

## BEIJING SEDECAL MEDICAL

Equipment Co., Ltd.  
104 Unit, No.1 Building, No.7  
Rongchang East Street,  
Beijing Economic-technological  
Development Area, Beijing, CHINA

+86 10 67806788  
[info@sedecal.cn](mailto:info@sedecal.cn)  
[service@sedecal.cn](mailto:service@sedecal.cn)



# SEDECAL

[www.sedecal.com](http://www.sedecal.com)

